

VESALIO

DRIVEN TO ADVANCE PATIENT CARE IN
VASCULAR OCCLUSION BY PROVIDING
PHYSICIANS SUPERIOR TECHNOLOGY DESIGNED
TO IMPROVE CLINICAL OUTCOMES

envast™

CHOOSE TO REMOVE


CATH-LAB STAFF
PRESENTATION



ENDOVASCULAR TREATMENT MODALITIES FOR STEMI


Standard endovascular treatment involves opening the artery with balloon & stent

- Outcomes remain poor in ~50% due to residual thrombus in the vessel



Endovascular aspiration continues to be debated



- CHEETAH, TOTAL, TASTE, TAPAS studies showed success in clot removal but also increased risk of AIS and variable outcomes



Studies show thrombus aspiration alone does not improve reperfusion or outcomes and bears higher potential for stroke

Neumann FJ, Sousa-Uva M, Ahlsson A, Alfonso F, Banning AP, Benedetto U, Byrne RA, Collet JP, Falk V, Head SJ, Juni P, Kastrati A, Koller A, Kristensen SD, Niebauer J, Richter DJ, Seferovic PM, Sibbing D, Stefanini GG, Windecker S, Yadav R, Zembala MO, Group ESCSD. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J 2019;40(2):87-165.


WHAT IS A STENT RETRIEVER



- A stent-retriever is a medical device
 - with a cylindrical mesh structure made of self-expanding nitinol
 - mounted on a wire
 - deployed within a micro-catheter
- Once at the site of the blood clot, the stent is released from within the catheter and self-expands within the thrombus
- This immediately pushes the clot against the wall of the artery, instantaneously reestablishing blood flow to the brain in most cases
- Standard of care for Ischemic Stroke



THE FIRST AND ONLY CE-APPROVED STENT-RETRIEVER
FOR CORONARY THROMBECTOMY



CHOOSE TO REMOVE

1. TREAT ALL LTB LESIONS

FROM SOFT CLOTS
THAT EASILY FRAGMENT
TO HARD, FIBRIN-RICH CLOTS
THAT CANNOT BE REMOVED

2. IMPROVE PROCEDURAL PERFORMANCE

A REAL SOLUTION FOR HIGH
CLOT BURDEN SITUATIONS

3. PROVIDE EASE OF USE

SYNERGISTIC WITH
ASPIRATION


Strive to achieve better patient outcomes

DROP ZONE™ TECHNOLOGY

A BALANCED DESIGN FOR SMOOTH TRACKING AND SAFE RETRIEVAL

DROP ZONES

entry points for large, organized thrombi



FLOW RESTORATION ZONE

radial force optimized for artery apposition

CLOSED DISTAL BASKET


clot retention inside structure

DESIGNED FOR
RAPID, HI-FLOW REPERFUSION

EFFICACY OUTCOMES – TIMI FLOW

enVast stent deployment was associated with immediate reperfusion in **85%** and TIMI-3 flow in **74%** of the patients with TIMI 0 after wire insertion


TIMI-3 increase from **31.7%** to **90%** after enVast ($p < .001$)






envast

NeVa Mechanical Thrombolysis System



envast Diameter **4.5 mm** Working Length **46 mm** Drop Zones **3** Wire Length **180 cm**



Contents: (1) envast 4.5 x 46 mm - flow restoration, 3 drop zones - 180 cm System
 Contenido: (1) Sistema envast de 4,5 x 46 mm - restauración del flujo, 3 zonas de captura - guía portadora de 180 cm
 Conteúdo: (1) Sistema envast de 4,5 x 46 mm - restabelecimento de fluxo, 3 zonas de captura - guia portadora de 180 cm

CATALOG NUMBER
REF EV-4546-F3RR LOT 051021X

Vesalio LLC
 105 North Pointe Drive
 Lake Forest, CA 92650
 +1-949-485-4708

ISO 13485
 MDSO Group
 Registered in
 9115-Farmington
 Connecticut

2023-05-10 Rx Only


REF EV-4546-F3RR	envast
LOT 051021X	4.5 mm / 46 mm
REF EV-4546-F3RR	envast
LOT 051021X	4.5 mm / 46 mm
REF EV-4546-F3RR	envast
LOT 051021X	4.5 mm / 46 mm
REF EV-4546-F3RR	envast
LOT 051021X	4.5 mm / 46 mm

Patents www.vesalio.com/patents LVLB-0104 REV B

envast Diameter **4.5 mm** Working Length **46 mm** Drop Zones **3** Wire Length **180 cm** REF EV-4546-F3RR LOT 051021X 2023-05-10



envast™ Diameter 4.5 mm Working Length 37 mm Drop Zones 2 Wire Length 180 cm



Contents: (1) envast 4.5 x 37 mm - Flow Restoration, 2 Drop Zones - 180 cm wire System
 Contenido: (1) Sistema envast de 4,5 x 37 mm – restauración del flujo, 2 zonas de captura – guía portadora de 180 cm
 Conteúdo: (1) Sistema envast de 4,5 x 37 mm – restabelecimento de fluxo, 2 zonas de captura – guia portadora de 180 cm

CATALOG NUMBER

REF EV-4537-F2RR LOT 050621X

REF EV-4537-F2RR
 LOT 050621X
 envast
 4.5 mm / 37 mm



 Vesalio LLC 105 North Pointe Drive Lake Forest, CA 92630 +615-206-7788	 Consult Instructions for Use Sterilized Using EO
	 Label Free
 MDSS GmbH Schiffgraben 41 30175 Hannover Germany	 Do Not Reuse
 Rx Only Prescription	 CE Mark 0297
 2023-05-06 Use By Date	 Keep Dry
	 Attention, Consult Accompanying Documents

REF EV-4537-F2RR
 LOT 050621X
 envast
 4.5 mm / 37 mm



REF EV-4537-F2RR
 LOT 050621X
 envast
 4.5 mm / 37 mm



REF EV-4537-F2RR
 LOT 050621X
 envast
 4.5 mm / 37 mm



Patents www.vesalio.com/patents

LV-LB-0102 REV B

envast™ Diameter 4.5 mm Working Length 37 mm Drop Zones 2 Wire Length 180 cm
 REF EV-4537-F2RR LOT 050621X 2023-05-06
 (01)00851279008712(17)230506(10)050621X

REF EV-4537-F2RR LOT 050621X 2023-05-06



PROCEDURE STEPS





CHOOSE TO REMOVE

enviast™



enVast™ CHOOSE TO REMOVE

- Choose an enVast size with labelled diameter that approximates the target vessel diameter

Product Name	Code	Maximal diameter	Working Length	Full Length	Drop Zones	Pusher Wire	Recommended Vessel Diameter (mm)	Min MC inner diameter
enVast 4.0 x 30	EV-4030-F2RR	4.0 mm	30 mm	39 mm	2	180 cm	≥ 2.0 & ≤ 3.5	.021"
								
enVast 4.5 x 37	EV-4537-F2RR	4.5 mm	37 mm	57 mm	2	180 cm	≥ 2.0 & ≤ 4.5	.021"
								
enVast 4.5 x 46	EV-4546-F3RR	4.5 mm	46 mm	66 mm	3	180 cm	≥ 2.0 & ≤ 4.5	.021"
								
enVast 6.0 x 35	EV-6035-F2RR	6.0 mm	35 mm	55 mm	2	180 cm	≥ 3.5 & ≤ 6.0	.027"
								

MICROCATHETER CONSIDERATIONS

- Choose a microcatheter size compatible with the enVast size chosen for the procedure

4.0 & 4.5 mm enVast sizes
are compatible with microcatheters with min ID of
0.021"

Via - 0.021"

Headway - 0.021"

TrevoPro - 0.021"

Phenom - 0.021"

Rebar 18 - 0.021"

Velocity - 0.025"

Marksman - 0.027"

Via - 0.027"

Phenom - 0.027"

6.0 mm enVast size
is compatible with microcatheters with min ID of 0.027"


Marksman - 0.027"

Via - 0.027"


Phenom - 0.027"

OTHER ACCESS CONSIDERATIONS


- Depending on the location of the lesion, you may choose to use:
 - A Guide Catheter + A Guide Catheter Extension: If clot is distal



- A Guide Catheter only: If clot is sufficiently proximal



ENVAST PREPARATION & FLUSHING




- Remove packaging hoop from the pouch and detach wire from the white rubber wire-holder
- Carefully remove delivery wire and introducer sheath out of the packaging hoop



- Flush enVast by inserting the distal end of the introducer sheath partially into the RHV connected to the microcatheter
- Tighten the RHV and verify that fluid exits the proximal end of the introducer sheath

LOADING ENVAST INTO THE MICROCATHETER



- Loosen the RHV and visually confirm that the tip of the sheath is seated deeply in the hub of the micro catheter



- Tighten the RHV around the introducer sheath to prevent back flow of blood, but not so tight as to damage the enVast device

NAVIGATING ENVAST TO THE OCCLUSION SITE

- Once the flexible portion of the pusher wire has entered the microcatheter shaft, loosen the RHV and remove the introducer sheath
- **DO NOT THROW AWAY THE INTRODUCER SHEATH:** Keep it on the sterile table as it may be needed for a second pass



NAVIGATING ENVAST TO THE OCCLUSION SITE

- Advance the enVast device in the microcatheter until the tip of the enVast basket aligns with the tip of the micro catheter under fluoro
- You can use the zebra markers on the enVast pusher wire to decide when to start visualizing the procedure under fluoroscopy
- If excessive resistance felt, discontinue delivery of enVast and identify the cause of resistance



ENVAST RETRIEVAL



- Use 3- 60ml VacLoc syringes.
- Check enVast position, take your time and start retrieval slowly
- Withdraw enVast and the GC-Extension simultaneously under continuous aspiration from the GC hub
- Aspirate GC, contrast run
- enVast can be used up to 3 times

WHAT DO YOU NEED TO HAVE AVAILABLE FOR STEMI THROMBECTOMY

ENVASt Thrombectomy Device



Remember to KEEP the INTRODUCER SHEATH!

2 x .014 Workhorse wires

Approved .021" or .027" ID Microcatheter

Guide Catheter and Guide Catheter Extension

3 -60 ml VacLoc Syringe
with 3-way stop cocks



enVast™ STEMI thrombectomy cheat sheet

ENVASt THROMBECTOMY DEVICE

Product Name	Code	Maximal diameter	Working Length	Full Length	Drop Zones	Pusher Wire	Recommended Vessel Diameter (mm)	Min MC inner diameter
enVast 4.0 x 30	EV-4030-F2RR	4.0 mm	30 mm	39 mm	2	180 cm	≥ 2.0 & ≤ 3.5	.021"
enVast 4.5 x 37	EV-4537-F2RR	4.5 mm	37 mm	57 mm	2	180 cm	≥ 2.0 & ≤ 4.5	.021"
enVast 4.5 x 46	EV-4546-F3RR	4.5 mm	46 mm	66 mm	3	180 cm	≥ 2.0 & ≤ 4.5	.021"
enVast 6.0 x 35	EV-6035-F2RR	6.0 mm	35 mm	55 mm	2	180 cm	≥ 3.5 & ≤ 6.0	.027"

LIST OF PRODUCTS REQUIRED

- enVast
- Approved .021 or .027 Microcatheter
- 3-60 ml VacLoc Syringe
- 3-way stop cocks
- Guide Catheter
- Guide Catheter Extension
- 2 x .014 Workhorse wires

MICROCATHETER CONSIDERATIONS

4.0 & 4.5 mm enVast sizes are compatible with MCs with min ID of 0.021" 6.0 mm enVast size is compatible with MCs with min ID of 0.027"

Rebar 18 - 0.021"	TrevoPro - 0.021"	Marksman - 0.027"
Headway - 0.021"	Via - 0.021"	Via - 0.027"
Phenom - 0.021"	Velocity - 0.025"	Phenom - 0.027"
Marksman - 0.027"	Via - 0.027"	
Phenom - 0.027"		

enVast™ STEMI thrombectomy cheat sheet


PROCEDURE STEPS

1. Cross occlusion site with .014" wire (GW-1) and go as distal as possible
2. Backload second .014 wire (GW-2) into .021" or .027" microcatheter (MC)
3. Advance MC and GW-2 to position MC tip past thrombus
4. Remove GW-2 from MC
5. Prep enVast device and advance into MC, drive to distal tip of MC
6. Unsheath enVast and remove the MC - Remember to KEEP the INTRODUCER SHEATH
7. If using a Guide Catheter Extension, advance over both enVast and GW-1
8. Attach 3 x 60 ml VacLoc Syringe with 3-way stop cocks to Guide Catheter (GC)
9. Check positioning and start slow retrieval
10. Withdraw enVast and GC-Extension simultaneously under continuous aspiration from the GC hub
11. Aspirate GC, perform fluoro run
12. If needed, repeat up to three enVast passes



en^vast™

QUESTIONS?



CHOOSE TO REMOVE