



DRIVEN TO ADVANCE PATIENT CARE
IN VASCULAR OCCLUSION BY
PROVIDING PHYSICIANS SUPERIOR TECHNOLOGY
DESIGNED TO IMPROVE CLINICAL OUTCOMES

NevaTM

Designed for 1st PASS SUCCESS with ALL Clot Types



LV-MKT-010-Rev B

WHY DEVELOP ANOTHER STENT-RETRIEVER

1

TREAT ALL OCCLUSIONS

FROM SOFT, FRIABLE CLOTS
THAT EASILY DISINTEGRATE
TO HARD, FIBRIN-RICH CLOTS
THAT ARE IMPENETRABLE

2

IMPROVE PROCEDURAL PERFORMANCE

1ST PASS SUCCESS
TIME TO RECANALIZATION
HIGHER TICI 2C/3 RATES

3

PROVIDE EASE OF USE

REAL TIME FEEDBACK
DURING RETRIEVAL
SYNERGISTIC WITH ALL
ACCESS PHILOSOPHIES

TO ACHIEVE BETTER PATIENT OUTCOMES



CONVENTIONAL STENT-RETRIEVERS



Work by **pinning** the clot to the artery wall and **dragging** it down

In most cases, clot **penetration** is **partial**

Hard clots simply slide outside the basket and **remain in place**

DESIGN


NevaTM

Designed for **1st PASS SUCCESS** with **ALL Clot Types**



NeVa™ DROP ZONE™ THE CLOT INSIDE

Uniquely designed to **CAPTURE**
ALL TYPES OF CLOT
INSIDE THE DEVICE
STRUCTURE




the **DROP ZONES™**
allow the capture of
large, organized thrombi
within the NeVa basket

NeVa™ A DESIGN TARGETING 1ST PASS SUCCESS

DROP ZONES™

2 or more Drop Zones offset at 90° work by acting as clot pockets: entry points to capture thrombi inside



BALANCED DESIGN

Optimized radial force
balanced with large
openings & closed ends

SMART MARKERS

2 per drop zone,
for real-time feedback
during retrieval

CLOSED DISTAL TIP


Clot gets inside,
clot stays inside!



THE NEXT FRONTIER IN THROMBECTOMY

the first and only
integrated Micro-Filtration technology

designed to
Maximize Clot Retention



The NeVa Net devices
incorporate a unique
braided net on the inside at
the distal end of the
expandable tip



5.5 x 37 mm

3 Drop Zones
VN-5537-03RR

Ideal for

**Proximal
occlusions**

Vessel diameters
3.5 – 5.5 mm

Recommended MC: 0.027"

4.5 x 29 mm

3 Drop Zones
VN-4529-03RR

Ideal for

**MCA
occlusions**

Vessel diameters
2.0 – 4.5 mm

Recommended MC > 0.021"

4.0 x 22 mm



2 Drop Zones
30020V-MS

Ideal for

**Distal M1, M2,
ACA, PCA occlusions**

Vessel diameters
2.0 – 3.5 mm

Recommended MC > 0.021"





5.5 x 37 mm

3 Drop Zones
VN-5537-03NR

Ideal for

**Proximal
occlusions**

Vessel diameters

3.5 – 5.5 mm

Recommended MC: 0.027"



4.0 x 30 mm

3 Drop Zones
VN-4030-03NR

Ideal for

**MCA
occlusions**

Vessel diameters

2.0 – 4.5 mm

Recommended MC > 0.027"

CLINICAL DATA

NevaTM

Designed for 1st PASS SUCCESS with ALL Clot Types



97% RECANALIZATION SUCCESS WITH 1.2 PASSES ACROSS ALL CLOT TYPES



Clot Type	Soft	Hard	Ultra Hard	All Clots
Clot morphology	Whole Blood "RED" Clot	Plasma Rich "WHITE" Clot	Clot modeled from ONYX 500	RED, WHITE and ONYX 500
N =	19	5	11	35
Length of clots - mm	10-40	6-12	4-12	4-40
1 st Pass TICI 3	84%	60%	55%	71%
Final TICI 3	89%	NR	82%	83%
Final TICI 2b/3	100%	100%	91%	97%
Average # of passes for final recanalization	1,05	1,00	1,63	1,23

CONSISTENT EFFECTIVENESS AT REMOVING ORGANIZED CLOTS

Data from Machi et al. Journal of Neuro-Int. Surgery, 2016¹


“All stent retrievers failed when interacting with large white thrombi (≥ 6mm)”

Solitaire*:	0/5	Trevo:	0/5
Embotrap*:	0/5	Eric:	0/5
Preset*:	0/5	Preset LT:	0/5
Catch*:	0/5	Separator 3D:	0/5
Revive*:	0/5	Mindframe:	0/5

Data from Machi P, et al., “Experimental evaluation of the NeVa™ thrombectomy device a novel stent retriever conceived to improve efficacy of organized clot removal”, Journal of Neuroradiology. 2018²

NeVa: 6/10 successful
complete removals
of white thrombi ≥ 6 mm

NeVa™ 1ST PASS RATES TRENDING HIGH



1. Abbasi M, Liu Y, Fitzgerald S, et al. Systematic review and meta-analysis of current rates of first pass effect by thrombectomy technique and associations with clinical outcomes. J Neurointerv Surg 2021;13:212-216

13 2. M. Ribo, et al., Journal of Neuroradiology May 2019

3. Bajrami et al. Interventional Neuroradiology Oct 2022



SUMMARY OF CLINICAL DATA

322 CASES WITH FIRST-LINE NEVA USE

	n	First Pass TICI 2C/3	First Pass TICI 2B/3	Final TICI 2C3	Final TICI 2B/3	Mean # of passes to achieve final TICI	Median # of passes to achieve final TICI (IQR)
Ribo et al.	30	14/30 47%	19/30 63%	19/30 63%	28/30 93%	1.7	
Akpinar et al.	118	53/118 45%	67/118 57%	86/118 73%	113/118 96%		1 (1-2)
Borggreffe et al.	29	14/29 48%	16/29 55%	21/29 72%	29/29 100%	2.4	
Bajrami et al.	145	78/145 54%	112/145 77%	126/145 87%	142/145 98%	1.8	
Total	322	159/322 49%	214/322 66%	252/322 78%	312/322 97%		

- M. Ribo, et al., Journal of Neuroradiology May 2019
- Akpinar, Cetin K., et al., Interventional Neuroradiology, July 2020
- Borggreffe et al. World Neurosurg. Dec. 2020
- Bajrami et al. Interventional Neuroradiology Oct 2022

1ST LINE TREATMENT ON “ALL” COMERS

30 PATIENTS -MULTI-CENTER

	2B/3	2C/3	
First Pass	19/30 → 63%	14/30 → 47%	Average # of passes for final recan → 1.7
Final Recanalization	28/30 → 93%	19/30 → 63%	

Patient Outcomes

- Mean NIHSS @ 24hr : 7
- 90 day mRS < 2 : 53%
- Zero NeVa related adverse events & sICH
- NeVa was effective with both balloon guide and local aspiration strategies
- In the 40 passes where the info was available:
70% clot incorporation into device basket



1ST LINE TREATMENT – ANTE & POSTERIOR

118 PATIENTS -MULTI-CENTER

	2B/3	2C/3	
First Pass	61/80 → 57%	46/80 → 45%	Median # of passes for final recan → 1 (IQR 1-2)
Final Recanalization	77/80 → 96%		

Favorable functional outcome (mRS ≤ 2):

- 53% in the “first-pass” subgroup
- 42.4% in the total patient population.

Procedure related complications:

- Symptomatic ICH: 3.3%
- Asymptomatic ICH: 13.6%
- Embolization into new territory: 1.7%
- Dissection that did not require stenting: 1.7 %

Flow Control strategies preferred:

- 92.4% of cases done with distal aspiration (Solumbra: Aspiration catheter + NeVa)
- Balloon Guide Catheter used only in 13.6% of cases



1ST LINE & RESCUE TREATMENT – M1

29 PATIENTS -SINGLE-CENTER

	2B/3	2C/3
First Pass	55%	48%
1-2 passes	79%	62%
Final Recanalization	100%	72%

Patient Outcomes:

- Median NIHSS scores decreased from 16 to 12 after treatment
- 90-day mRS < 2 : 31%

Safety Data:

- 1 asymptomatic carotico-cavernous fistula
- 1 asymptomatic M2 dissection
- 1 sICH
- Procedure-related vasospasm rate: 48% (no negative impact observed on outcomes)

1ST LINE TREATMENT – ANTERIOR

145 PATIENTS -SINGLE-CENTER

	2B/3	2C/3	
First Pass	112/145 → 77%	78/145 → 54%	Average # of passes for final recan → 1.8
Final Recanalization	142/145 → 98%	126/145 → 87%	

Patient Outcomes:

- 24 hour mean NIHSS:
 - 6 (IQR 0-33)
 - average decrease of 10 points compared to initial presentation
- 30-day mRS 0-2: 89 (60.7%)
- 90-day mRS 0-2: 97 (66.9%)

Safety Data:

- No device-related adverse events
- No significant vasospasm
- HT-1: 29.7% (n=43), HT-2: 7.6% (n=11), PH-1: 2.8% (n=4)
- no perforation or extravasation observed during any procedure
- asymptomatic SAH identified on follow-up imaging: 4.8%(n=7)

EVALUATION

NevaTM

Designed for 1st PASS SUCCESS with ALL Clot Types



Proposal:

Retrospective analysis of 10 consecutive incoming AIS patients

No particular patient exclusion criteria, hospital protocol to be followed, but recommend to start with standard cases to gain familiarity with NeVa

NeVa tips & tricks training before use

Expectations:

Use NeVa as first line treatment

At least 3 attempts to achieve TICl 2b/3 before trying an alternative device

A simple form to fill for each case

THANK YOU

NevaTM

Designed for 1st PASS SUCCESS with ALL Clot Types



SUPPLEMENTARY INFORMATION

RADIAL FORCE

NEVA MARKERS


EXAMPLE CASES &
CLOTS BY NEVA

NeviaTM

Designed for 1st PASS SUCCESS with ALL Clot Types




Neva™ DROP ZONE™ THE CLOT INSIDE




OPTIMIZED RADIAL FORCE BALANCED WITH LARGE OPENINGS & CLOSED ENDS

Expansive Radial Force



Compressive Radial Force




— NeVa 5.5 x 37 mm — NeVa 4.0 x 22 mm — NeVa 4.5 x 29 mm — Solitaire 6x40 mm

— NeVa 5.5 x 37 mm — NeVa 4.0 x 22 mm — NeVa 4.5 x 29 mm — Solitaire 6x40 mm

THE LARGEST NEVA SIZE (5.5) RADIAL FORCE COMPARED TO SIMILAR SIZED STENT RETRIEVER

Compressive Radial Force Measurements





A PROXIMAL OCCLUSION: ONE AND DONE

Right ICA Tip Occlusion, 1st Pass Success

NeVa 4.5 x 37 mm

Prof Geyik, Aydin University, Istanbul, TURKEY

[READ CASE STUDY >](#)



1ST PASS IN BASILAR WAKE UP STROKE

Basilar Occlusion, 1st Pass Success

NeVa 4.5 x 29 mm

Dr Sirvinskas, Republic University, Vilnius, LITHUANIA

[READ CASE STUDY >](#)



1ST PASS SUCCESS WITH 3 DROP ZONES

Left M1 Occlusion, first pass success

NeVa 4.5 x 29 mm

Dr Maurer, University Hospital, Augsburg, GERMANY

[READ CASE STUDY >](#)



1ST PASS SUCCESS AFTER CAROTID BLOWOUT REPAIR

Left M2 Occlusion, first pass success through the carotid stent graft

NeVa 4.0 x 22 mm

Prof Kizilkilic, Dr Korkmaz, Cerrahpasa University, Istanbul, TURKEY

[READ CASE STUDY >](#)



1ST PASS IN STROKE WITH UNKNOWN ONSET

Right M1 Occlusion, 1st Pass Success

NeVa 4.0 x 30 mm

Dr. Kalousek, Sisters Charity Hospital, Zagreb, Croatia

[READ CASE STUDY >](#)



NEVA IN TANDEM STROKE

Tandem Occlusion, two single-pass retrievals, case from LINNC MASTERCLASS

NeVa 4.0 x 30 mm

Prof Spelle, Prof Moret, Dr Mihalea, Neuri Bicetre, Paris, FRANCE

[WATCH CASE >](#)



IMPACT OF 1ST PASS SUCCESS IN EARLY ONSET STROKE

Left M1 Occlusion, first pass success

NeVa 4.0 x 30 mm

Prof Mayer, University Hospital, Jena, GERMANY

[READ CASE STUDY >](#)



WAKE UP STROKE 1ST PASS SUCCESS

Left M1 Occlusion, 1st Pass Success

NeVa 4.0 x 30 mm

Prof Geyik, Aydin University Hospital, Istanbul, TURKEY

[READ CASE STUDY >](#)



SINGLE NEVA RESCUES KISSING RETRIEVERS

Carotid T Occlusion, 1st Pass Success after 2 failed attempts with the kissing-stents technique

NeVa 6.0 x 44 mm

Dr Tomasello, Vall d'Hebron, Barcelona, SPAIN

[READ CASE STUDY >](#)



NEVA TO THE RESCUE

Left M2 Occlusion, single pass rescue after failure of 2 different devices

NeVa 4.0 x 22 mm

Prof Geyik, Aydin University, Istanbul, TURKEY

[READ CASE STUDY >](#)



NEW! NEVA SAVES THE DAY AFTER A 5-PASS ORDEAL

Left M1 Occlusion Success

NeVa™ 4.5 x 37 mm

Bucharest University Emergency Hospital Stroke Team

[READ CASE STUDY >](#)



NEW! NEVA IN AN I-TYPE ICA WITH MANY SURPRISES

Live case transmission from iCureStroke 2022

NeVa™ 4.5 x 37 mm

Prof Geyik & Dr Bajrami, Aydin University Hospital, Istanbul, TURKEY

[WATCH CASE >](#)

Ne^va™ DROP ZONE™ THE CLOT INSIDE

