



CHANGING OUTCOMES CHANGING LIVES

Designed for 1st PASS SUCCESS with ALL Clot Types



Neva™


VESALIO CHANGING OUTCOMES CHANGING LIVES

2017 Founded by physicians treating stroke

- NeVa design freeze and establishment of Vesalio in 2017
- 34 Issued & 50+ Filed Patents

2018 Set on resolving vascular occlusions

- Commercial launch of NeVa in Europe



2022 Improving, perfecting, diversifying portfolio

- U.S. FDA Approval for Vasospasm (NeVa VS)
- CE marking of NeVa NET - the 1st SR device with integrated distal filter in thrombectomy
- CE marking of enVast - the 1st SR-type device approved in STEMI

2023 Commercial Expansion and Success

- International commercialization in over 50 countries, expanding into new global regions
- U.S. commercialization with NeVa VS
- 10000th device milestone

Vesalio is advancing the care of patients suffering from vascular occlusion by providing physicians superior technology designed to improve clinical outcomes

WHY DEVELOP ANOTHER STENT-RETRIEVER

1

TREAT ALL OCCLUSIONS

FROM SOFT, FRIABLE CLOTS
THAT EASILY DISINTEGRATE
TO HARD, FIBRIN-RICH CLOTS
THAT ARE IMPENETRABLE

2

IMPROVE PROCEDURAL PERFORMANCE

1ST PASS SUCCESS
TIME TO RECANALIZATION
HIGHER TICI 2C/3 RATES

3

PROVIDE EASE OF USE

REAL TIME FEEDBACK
DURING RETRIEVAL
SYNERGISTIC WITH ALL
ACCESS PHILOSOPHIES

TO ACHIEVE BETTER PATIENT OUTCOMES

CONVENTIONAL STENT-RETRIEVERS



Work by **pinning** the clot to the artery wall and **dragging** it down

In most cases, clot **penetration** is **partial**

Hard clots simply slide outside the basket and **remain** in place

DESIGN


NevaTM

Designed for 1st PASS SUCCESS with ALL Clot Types



NeVa™ DROP ZONE™ THE CLOT INSIDE

Uniquely designed to
CAPTURE ALL TYPES OF CLOT
INSIDE THE DEVICE STRUCTURE



the DROP ZONES™
allow the capture of
large, organized thrombi
within the NeVa basket

DROP ZONES™

2 or more Drop Zones offset at 90° work by acting as clot pockets: entry points to capture thrombi inside



BALANCED DESIGN


Optimized radial force balanced with large openings & closed ends

DROP ZONE MARKERS

2 per Drop Zone, for real-time feedback during retrieval


CLOSED DISTAL TIP

Clot gets inside, clot stays inside!




OPTIMIZED RADIAL FORCE BALANCED WITH LARGE OPENINGS & CLOSED ENDS

Expansive Radial Force



Compressive Radial Force



— NeVa 5.5 x 37 mm — NeVa 4.0 x 22 mm — NeVa 4.5 x 29 mm — Solitaire 6x40 mm


— NeVa 5.5 x 37 mm — NeVa 4.0 x 22 mm — NeVa 4.5 x 29 mm — Solitaire 6x40 mm





NEVA 5.5 COMPRESSIVE RADIAL FORCE COMPARED TO SOLITAIRE 6.0

Compressive Radial Force Measurements





5.5 x 37 mm

3 Drop Zones
VN-5537-03RR

Ideal for

**Proximal
occlusions**

Vessel diameters

3.5 – 5.5 mm

Recommended MC: 0.027"

4.5 x 29 mm

3 Drop Zones
VN-4529-03RR

Ideal for

**MCA
occlusions**

Vessel diameters

2.0 – 4.5 mm

Recommended MC > 0.021"

4.0 x 22 mm

2 Drop Zones
30020V-MS

Ideal for

**Distal M1, M2,
ACA, PCA occlusions**

Vessel diameters

2.0 – 3.5 mm

Recommended MC > 0.021"


CLINICAL DATA

NevaTM

Designed for 1st PASS SUCCESS with ALL Clot Types



97% RECANALIZATION SUCCESS WITH 1.2 PASSES ACROSS ALL CLOT TYPES



Clot Type	Soft	Hard	Ultra Hard	All Clots
Clot morphology	Whole Blood "RED" Clot	Plasma Rich "WHITE" Clot	Clot modeled from ONYX 500	RED, WHITE and ONYX 500
N =	19	5	11	35
Length of clots - mm	10-40	6-12	4-12	4-40
1 st Pass TICI 3	84%	60%	55%	71%
Final TICI 3	89%	NR	82%	83%
Final TICI 2b/3	100%	100%	91%	97%
Average # of passes for final recanalization	1,05	1,00	1,63	1,23

CONSISTENT EFFECTIVENESS AT REMOVING ORGANIZED CLOTS

Data from Machi et al. Journal of Neuro-Int. Surgery, 2016¹

“All stent retrievers failed when interacting with large white thrombi (≥ 6mm)”

Solitaire*:	0/5	Trevo:	0/5
Embotrap*:	0/5	Eric:	0/5
Preset*:	0/5	Preset LT:	0/5
Catch*:	0/5	Separator 3D:	0/5
Revive*:	0/5	Mindframe:	0/5

Data from Machi P, et al., “Experimental evaluation of the NeVa™ thrombectomy device a novel stent retriever conceived to improve efficacy of organized clot removal”, Journal of Neuroradiology. 2018²

**NeVa: 6/10 successful
complete removals
of white thrombi ≥ 6 mm**



1ST PASS RATES TRENDING HIGH

CLINICAL RESULTS

META-ANALYSIS¹
published in 2021


16870 patients, 67 clinical series

NeVa CLEAR study
Presented in 2023²

107 patients (mITT), 25 centers

NeVa largest published
patient series in 2022³

145 patients, single center



1. Abbasi M, Liu Y, Fitzgerald S, et al. Systematic review and meta-analysis of current rates of first pass effect by thrombectomy technique and associations with clinical outcomes. J Neurointerv Surg. 2021;13:212-216.
2. Yoo AJ, Geyik S, Froehler MT, et al. Primary results from the CLEAR study of a novel stent retriever with drop zone technology. Journal of NeuroInterventional Surgery. Published Online First: 02 December 2023. doi: 10.1136/jnis-2023-020960
3. Bajrami A, Ertugrul O, Senadim S, Erdem E, Baltacioglu F, Geyik S. First pass results of mechanical thrombectomy with two-drop zone NeVaTM device. Interv Neuroradiol. 2022 Oct 30;15910199221135309. doi: 10.1177/15910199221135309. PMID: 36314456.



PRIMARY RESULTS FROM THE CLEAR STUDY OF A NOVEL STENT RETRIEVER WITH DROP ZONE TECHNOLOGY

Albert J Yoo, Texas Stroke Inst, Plano, TX

Christoph Maurer, Augsburg, Augsburg, Germany

Raul G Nogueira, Univ of Pittsburgh, Pittsburgh, PA

Laurent Spelle, Bicetre Univ Hosp, Paris, France

Audrius Sirvinskas, Vilnius, Vilnius, Lithuania

Jordi Blasco, Hosp Clinic Barcelona, Barcelona, Spain

Adnan H Siddiqui, Univ at Buffalo, Buffalo, NY

Diogo C Haussen, Emory Univ, Atlanta, GA

Arturo Consoli, Hosp Foch, Paris, France

Alberto Maud, Texas Tech Univ, El Paso, TX

Maher Sahnoun, CHU Reims, Reims, France

Serdar Geyik, Istanbul Aydin Univ, Istanbul, Turkey

Tareq Kass-Hout, Univ of Chicago, Chicago, IL

Ricardo Hanel, Baptist, Jacksonville, FL

Demetrius Lopes, Advocate Aurora Health, Park Ridge, IL

Eugene Lin, Mercy St. Vincent, Toledo, OH

Muhammad Asif Taqi, Vasc. Neurology of S.Cal, Thousand Oaks, CA

Timothy Ryan Miller, Univ of Maryland, Baltimore, MD

Keith Woodward, Fort Sanders, Fort Sanders, TN

M. Imran Chaudry, Greenville, Mt. Pleasant, SC

Joshua Bentley, Southeast Health, Dothan, AL

Jens Fiehler, Uni Medical Ctr Hamburg-Eppendorf, Hamburg, Germany

Michael Froehler, Vanderbilt, Nashville, TN

Osama Zaidat, Mercy Health St. Vincent Hosp, Toledo, OH

Laurent Pierot, CHU Reims, Reims, France

Ameer E Hassan, Valley Baptist Medical Ctr, Harlingen, TX

Marc Ribo, Val d'Hebron, Barcelona, Spain

Aamir Badruddin, Community Healthcare, Munster, IN

Shazam Muhammad Hussain, Cleveland Clinic, Cleveland, OH

Christoph Groden, Univsmedizin Mannheim, Mannheim, Germany

Christian Ramsey, Baptist Health, Lexington, KY

Arsida Bajrami, Istanbul Aydin Univ, Istanbul, Turkey

Rishi Gupta, Wellstar Health System, Marietta, GA

NevaTM

Journal of Interventional Surgery

Published Online First: 02 December 2023. doi: 10.1136/jnis-2023-020960

Yoo AJ, Geyik S, Froehler MT, et al Primary results from the CLEAR study of a novel stent retriever with drop zone technology. Journal of NeuroInterventional Surgery. Published Online First: 02 December 2023. doi: 10.1136/jnis-2023-020960



OBJECTIVES

1. Demonstrate the efficacy and safety of the NeVa device for revascularization of LVOs
2. Assess the comparative performance of the NeVa device versus existing stent retriever devices

STUDY DESIGN

Prospective, multi-center, open label, single arm, FDA-regulated IDE study

POPULATION:

- AIS patients treatable within 8 hours of onset
- 18-85 yrs
- NIHSS 8-25
- ASPECTS 6-10
- intracranial LVO including ICA, M1/M2 MCA, and vertebrobasilar arteries
- IV-tPA ≤ 3 h of TLKW



STUDY ENDPOINTS

PRIMARY EFFICACY ENDPOINT

- Rate of successful reperfusion (eTICI 2b-3) within 3 NeVa passes without rescue
- Primary non-inferiority analysis (mITT population)
- Performance goal: 72% with non-inferiority margin of -10%

PRIMARY SAFETY ENDPOINT

- All-cause mortality at 90 days and/or symptomatic ICH (sICH) at 24 hours post procedure (ECASS III)

SECONDARY ENDPOINTS: (HIERARCHICAL TESTING)

- First pass eTICI 2b-3 with NeVa
- eTICI 2b-3 after all NeVa passes (no rescue)
- eTICI 2c-3 after all NeVa passes (no rescue)
- 90-day good outcome (mRS 0-2)

METHODOLOGY

- mITT population was used for the analysis of efficacy endpoints:
 - All subjects that met eligibility criteria and used the 4.0x22mm, 4.5x29mm or the 5.5x37mm NeVa devices
 - N=107
- ITT population was used for the analysis of safety endpoints
 - All subjects that got treated with the NeVa devices, including those that did not meet eligibility criteria
 - N=139
- Performance Goal (PG): defined using ARISE 2 and TIGER trials

SUMMARY OF RESULTS

PARAMETER	ITT/ SAFETY (n = 139)	mITT (n = 107)
Age, mean (SD)	66.7 (12.8)	65.1 (13.2)
Baseline NIHSS, median (IQR)	16 (12-20)	16 (12-20)
Baseline CT ASPECTS, median (IQR)	9 (8-10)	9 (8-10)
IV tPA administration, n (%)	71 (51.1%)	59 (55.1%)

PROCEDURAL CHARACTERISTICS

Last known well to arterial puncture min; median (IQR)	202 (138-294)	181 (131-252)
Arterial puncture to first device pass min; median (IQR)	18 (11-25)	18 (11-25)
Procedure duration min; median (IQR)	32 (20-51)	35 (22-52)


PRIMARY OCCLUSIVE LOCATION

Intracranial ICA, n (%)	14 (10.1%)	10 (9.3%)
MCA-M1, n (%)	86 (61.9%)	66 (61.7%)
MCA-M2, n (%)	37 (26.6%)	30 (28.0%)
Basilar, n (%)	1 (0.7%)	1 (0.9%)
PCA, n (%)	1 (0.7%)	0 (0.0%)


USE OF ANCILLIARY DEVICES FOR FLOW CONTROL

■ IC Only ■ BGC Only ■ BGC + IC

ITT




mITT



eTICI 2B-3 WITHIN 3 PASSES

Primary efficacy endpoint (mITT, n=107)



NeVa successfully demonstrated non-inferiority to the performance goal based on predicate devices (prespecified primary efficacy analysis)


NeVa further demonstrated post hoc superiority to the predicate performance goal

FIRST-PASS REPERFUSION


Key secondary efficacy endpoint (mITT, n=107)

NeVa achieved a superior rate of first pass successful reperfusion compared to predicate devices - prespecified secondary endpoint hierarchical testing

SUCCESSFUL RECANALIZATION % eTICI 2B-3



EXCELLENT RECANALIZATION % eTICI 2C-3



FIRST-PASS REPERFUSION

Key secondary efficacy endpoint (mITT, n=107)

NeVa achieved a superior rate of first pass successful reperfusion compared to predicate devices - prespecified secondary endpoint hierarchical testing

SUCCESSFUL RECANALIZATION
% eTICI 2B-3

≈ 7.5
10

EXCELLENT RECANALIZATION
% eTICI 2C-3

≈ 5
10

OTHER

ADDITIONAL ANGIOGRAPHIC EFFICACY ENDPOINTS

Successful Reperfusion after all passes (eTICI 2B-3) – n, %
Excellent Reperfusion after all passes (eTICI 2C-3) – n, %

mITT (n = 107)	PREDICATE DEVICE VALUES	P VALUE
102 (95.3%)	95.7%	0.85
76 (71.0%)	71.8%	0.86


PROCEDURAL STATISTICS

Use of Rescue Therapy – n, %
Time from arterial puncture to 1st device pass (min), median (IQR)
Procedure Duration (min), median (IQR)
Number of Passes, median (IQR)

4 (3,7%)
18 (11 – 25)
32 (20 – 51)
1 (1 – 2)

COMPOSITE OF 90-DAY ALL-CAUSE MORTALITY AND/OR 24-H SICH

Primary safety endpoint (ITT, n=139)



OTHER SAFETY ENDPOINTS

Percent of subjects	ITT	95% Confidence Interval
deceased at Day 90	9.4%	5.5% - 15.3%
experiencing sICH at 24 hrs post-procedure	5.0%	2.5% - 10.0%
with > 4-point increase in NIHSS score at 24 hrs post-procedure	5.0%	2.5% - 10.0%
with > 4-point increase in NIHSS score at day 5-10/ Discharge	2.9%	1.1% - 7.2%
experiencing PRSAEs (procedure-related serious adverse events)	4.3%	
experiencing DRSAEs (device-related serious adverse events)	3.6%	
		<u>NeVa-related SAEs:</u> Vasospasm (2), aSAH (1), SAH (2), vessel perforation (1)


CLINICAL OUTCOMES

90-day outcomes with NeVa were superior to predicate studies in exploratory analysis

90-DAY MRS ≤ 2

(mITT, n=106)

65.1%
vs. PG: 54.7%, P = 0.03



OTHER CLINICAL OUTCOMES

(mITT, n=107)

NIHSS at 24-hours, median (IQR) 4 (1-8)


NIHSS change from baseline to 24-hrs, median (IQR) -10 (-15 to -5)

NIHSS at 5-10 days or discharge, median (IQR) 2 (0-5)

NIHSS change from baseline to 5-10 days or discharge, median (IQR) -11 (-17 to -8)

CONCLUSION

- The CLEAR study demonstrates that the NeVa stent retrievers are effective and safe for the revascularization of LVOs in AIS
 - Substantial equivalence to predicate devices (eTICI 2b-3) in primary non-inferiority analysis
 - Statistically superior to predicate devices (eTICI 2b-3) in post hoc analysis
- Statistically superior to predicate devices for first pass successful reperfusion (eTICI 2b-3)
- Clinical outcomes at 90 days were superior to predicate studies in exploratory analysis



EVALUATION

NevaTM

Designed for 1st PASS SUCCESS with ALL Clot Types



Proposal:

Retrospective analysis of 10 consecutive incoming AIS patients

No particular patient exclusion criteria, hospital protocol to be followed, but recommend to start with standard cases to gain familiarity with NeVa

NeVa tips & tricks training before use

Expectations:

Use NeVa as first line treatment

At least 3 attempts to achieve TICl 2b/3 before trying an alternative device

A simple form to fill for each case



A PROXIMAL OCCLUSION: ONE AND DONE

Right ICA Tip Occlusion, 1st Pass Success

NeVa 4.5 x 37 mm

Prof Geyik, Aydin University, Istanbul, TURKEY

[READ CASE STUDY >](#)



1ST PASS IN STROKE WITH UNKNOWN ONSET

Right M1 Occlusion, 1st Pass Success

NeVa 4.0 x 30 mm

Dr. Kalousek, Sisters Charity Hospital, Zagreb, Croatia

[READ CASE STUDY >](#)



SINGLE NEVA RESCUES KISSING RETRIEVERS

Carotid T Occlusion, 1st Pass Success after 2 failed attempts with the kissing-stents technique

NeVa 6.0 x 44 mm

Dr Tomasello, Vall d'Hebron, Barcelona, SPAIN

[READ CASE STUDY >](#)



1ST PASS IN BASILAR WAKE UP STROKE

Basilar Occlusion, 1st Pass Success

NeVa 4.5 x 29 mm

Dr Sirvinskas, Republic University, Vilnius, LITHUANIA

[READ CASE STUDY >](#)



NEVA IN TANDEM STROKE

Tandem Occlusion, two single-pass retrievals, case from LINNC MASTERCLASS

NeVa 4.0 x 30 mm

Prof Spelle, Prof Moret, Dr Mihalea, Neuri Bicetre, Paris, FRANCE

[WATCH CASE >](#)



NEVA TO THE RESCUE

Left M2 Occlusion, single pass rescue after failure of 2 different devices

NeVa 4.0 x 22 mm

Prof Geyik, Aydin University, Istanbul, TURKEY

[READ CASE STUDY >](#)



1ST PASS SUCCESS WITH 3 DROP ZONES

Left M1 Occlusion, first pass success

NeVa 4.5 x 29 mm

Dr Maurer, University Hospital, Augsburg, GERMANY

[READ CASE STUDY >](#)




IMPACT OF 1ST PASS SUCCESS IN EARLY ONSET STROKE

Left M1 Occlusion, first pass success

NeVa 4.0 x 30 mm

Prof Mayer, University Hospital, Jena, GERMANY

[READ CASE STUDY >](#)



NEW! NEVA SAVES THE DAY AFTER A 5-PASS ORDEAL

Left M1 Occlusion Success

NeVa™ 4.5 x 37 mm

Bucharest University Emergency Hospital Stroke Team

[READ CASE STUDY >](#)



1ST PASS SUCCESS AFTER CAROTID BLOWOUT REPAIR

Left M2 Occlusion, first pass success through the carotid stent graft

NeVa 4.0 x 22 mm

Prof Kizilkilic, Dr Korkmaz, Cerrahpasa University, Istanbul, TURKEY

[READ CASE STUDY >](#)



WAKE UP STROKE 1ST PASS SUCCESS

Left M1 Occlusion, 1st Pass Success

NeVa 4.0 x 30 mm

Prof Geyik, Aydin University Hospital, Istanbul, TURKEY

[READ CASE STUDY >](#)



NEW! NEVA IN AN I-TYPE ICA WITH MANY SURPRISES

Live case transmission from iCureStroke 2022

NeVa™ 4.5 x 37 mm

Prof Geyik & Dr Bajrami, Aydin University Hospital, Istanbul, TURKEY

[WATCH CASE >](#)



THANK YOU

Designed for 1st PASS SUCCESS with ALL Clot Types

NevaTM

