



A FIRST PASS TICI 3 IN THE RIGHT POSTERIOR CEREBRAL ARTERY

NeVa™ 4.0 x 22 mm

Dr. Lucian Gramada, Dr. Huda Elmoula
Universitätsklinikum Witten/ Herdecke, Germany

PRESENTATION

85-year-old male patient
presented with a NIHSS of 2
within the mechanical
thrombectomy eligibility time
window.

No IV-tPA was administered.





Dr. Lucian Gramada
Universitätsklinikum
Witten/Herdecke
Witten, Germany



Dr. Huda Elmoula
Universitätsklinikum
Witten/Herdecke
Witten, Germany

PRE-OP IMAGING


Imaging confirmed occlusion of right PCA and a slight stenosis.



THE PROCEDURE

- The lesion was accessed via the right femoral artery.
- Dr. Gramada and Dr. Elmoula used a NeVa 4.0 x 22 mm and achieved a TICl 3 level recanalization after a single pass using a combined technique with NeVa and a 0.062" ID aspiration catheter.
- The mild stenosis didn't need further treatment.
- Recanalization was achieved within 36 minutes from groin puncture.

ANGIOGRAPHIC IMAGING BEFORE, DURING, AND AFTER THE FIRST THROMBECTOMY PASS




ANGIOGRAPHIC IMAGING AT THE BEGINNING OF THE PROCEDURE

Confirmed the occlusion of the right PCA branch.



ANGIOGRAPHIC IMAGING DURING THE PROCEDURE

Showing NeVa 4.0 x 22 mm deployed at the occlusion site.



ANGIOGRAPHIC IMAGING AT THE END OF THE PROCEDURE

Showing TICl 3 recanalization after the first NeVa pass.



ANGIOGRAPHIC IMAGING AT THE END OF THE PROCEDURE

Complete recanalization of the right-PCA was possible in a single pass.

CLINICAL OUTCOME

24-hour NIHSS: 0
24-hour mRS: 0

 DR HUDA ELMOULA

“Nice try stroke, but I’m in the game.”