

Endovascular TODAY

LEVELING UP IN LIMB PRESERVATION

Multifaceted
approaches to
delivering truly
meaningful and
better outcomes
in CLTI.

MODERN TRIAL
ENDPOINTS

CLINICAL
STUDY NEEDS

AWARENESS
INITIATIVES

CENTERS OF
EXCELLENCE

ADVANCED
TECHNIQUES

pVasc™ Thrombectomy System: A Simple and Elegant Solution for Thromboembolic Disease

pVasc offers the lowest-profile 0.035-inch catheter-compatible delivery to remove diverse clot morphologies from the peripheral arterial vasculature.

By Halim Yamine, MD

Acute limb ischemia (ALI) is a severe limb-threatening condition. The most common causes include a local thrombotic event from a preexisting plaque or an embolus from a proximal source such as the heart. ALI requires prompt diagnosis and treatment due to a high risk of limb loss (15%-20%) and mortality (15%-20%).^{1,2}

Treatment of ALI includes percutaneous interventions such as catheter-directed thrombolysis, thromboaspiration, pharmacomechanical thrombectomy, and mechanical thrombectomy, in addition to traditional surgical interventions. Given the urgent nature of ALI, there are important considerations in selecting the appropriate treatment modality, including time to

"pVasc has become an important tool in my practice. It is very easy to use and allows me to treat multiple vascular beds, including distal tibial arteries, with a single device and minimal blood loss. It is very well designed and does well in retrieving both soft thrombus and more organized clot."

– Dr. Yamine




Figure 1. pVasc key features and benefits.




Figure 2. Diverse clot morphologies captured by pVasc.

reperfusion, location of the lesion, and the frequent need to treat small tibial vessels.

The pVasc Thrombectomy System (Vesalio, Inc.) offers a simple and effective percutaneous mechanical solution for removing thromboembolus from vessels ranging 2 to 6 mm in diameter (Figure 1). Comprised of a self-expanding nitinol structure mounted on a 200-cm pusher wire, pVasc's low-profile design allows delivery through a 0.035-inch wire-compatible catheter,* providing exceptional navigation to distal arteries.

pVasc features a streamlined, intuitive design that enables easy handling to minimize device manipulation for efficient delivery and clot retrieval with minimal blood loss. It is engineered with proprietary Drop Zone™ technology (Vesalio, Inc.) to trap and retain thromboembolus inside the device's lumen for secure clot removal. Combined with a closed distal basket, pVasc is designed to address diverse clot morphologies and organized emboli while minimizing the risk of clot migration (Figure 2). There is no capital equipment involved.

CASE REPORT

Case Presentation

A woman in her late 70s presented with acute left lower extremity ischemia due to popliteal artery occlusion from embolization.

After contralateral access was achieved, a 7-F, 90-cm sheath was introduced into the proximal left popliteal artery. An initial angiogram revealed complete occlusion of the popliteal artery as well as the proximal tibial arteries with reconstitution (Figure 3).

Procedural Overview

A 6-mm pVasc device was delivered through a 0.035-inch CXI catheter (Cook Medical) into each of the three tibial vessels for clot removal.

In each pass, the CXI catheter was introduced into the tibial vessel over a wire and placed distal to the clot such that each Drop Zone was fully utilized to maximize clot capture. Next, pVasc was advanced into the catheter until its tip reached the distal end of the catheter without exiting it. pVasc was deployed by retracting the catheter until the device was fully exposed. With pVasc deployed and firmly in place, the catheter was then completely removed before withdrawing the device slowly into the sheath. The sheath valve was removed before pVasc exited the sheath. After pVasc was completely removed, the sheath valve was then flushed and reinserted to perform follow-up angiography. Finally, the device was cleaned on the back table. These steps were repeated in the next vessel.

We began with pVasc in the peroneal artery (Figure 4), followed by the anterior tibial artery (ATA) (Figure 5), and finally the posterior tibial artery (PTA) (Figure 6). Although not used in this case, a 0.014-inch wire can be maintained alongside the 0.035-inch CXI catheter as a buddy wire to facilitate reaccess.

pVasc easily navigated into each vessel and successfully removed clot with each pass. The device took the curve of the ATA very easily and withdrew smoothly without any issues.

Case Conclusion

Complete flow restoration was achieved (Figure 7) with three-vessel runoff all the way to the foot in

pVasc THROMBECTOMY SYSTEM

Sponsored by Vesalio, Inc.




Figure 3. Initial angiogram.



Figure 4. pVasc deployed in the peroneal artery.



Figure 5. pVasc deployed in the ATA.



Figure 6. pVasc delivered through 0.035-inch catheter into the PTA.




Figure 7. Final angiogram.




Figure 8. Clot secured on pVasc (A) and samples of captured clot (B).

< 60-minute total procedure time from skin to skin. Pulses were present in the dorsalis pedis artery and PTA post-intervention. Total blood loss was about 20 mL.


CONCLUSION

Since its release in the fall of 2024, the pVasc Thrombectomy System has been used successfully to retrieve thromboembolus of varying composition with its proprietary Drop Zone™ technology. Early feedback highlights pVasc's low-profile design, which provides exceptional navigability through a 0.035-inch catheter for traversing stenotic vessels and working in distal arteries. Along with the device's easy handling

and simple setup, pVasc offers a powerful, streamlined thrombectomy solution to answer the challenges in ALI treatment. ■

**Minimum delivery catheter inner diameter for 4-mm pVasc is 0.021 inches and for 6-mm pVasc is 0.027 inches.*

1. Siah MC, Shih M. Evaluation and management of acute limb ischemia. *Endovasc Today*. 2021;20:63-66. <https://evtoday.com/articles/2021-jan/evaluation-and-management-of-acute-limb-ischemia>
2. Obara H, Matsubara K, Kitagawa Y. Acute limb ischemia. *Ann Vasc Dis*. 2018;11:443-448. doi: 10.3400/avd.ra.18-00074



Halim Yammine, MD

Clinical Assistant Professor of Vascular Surgery
Wake Forest School of Medicine
Associate Program Director
Vascular Surgery Fellowship
Sanger Heart and Vascular Institute
Atrium Health
Charlotte, North Carolina
halim.yammine@atriumhealth.org
Disclosures: Consultant to Vesalio, Inc., Penumbra, and Teleflex; receives an educational grant from Medtronic.

