

DO THE DROP ZONE

Discover NeVa and NeVa NET

NeVaTM


Designed for 1st PASS SUCCESS with ALL Clot Types



LV-MKT-010 REV E

VESALIO CHANGING OUTCOMES CHANGING LIVES

Vesalio is advancing patient care in vascular occlusion by innovating in thrombectomy technologies to improve clinical outcomes



Founded by
physicians treating
vascular occlusion

Backed by world-class
investors

Present in more than
60 countries
worldwide

TECHNICAL CHALLENGES OF THROMBECTOMY

NeVa was designed to address the technical challenges of thrombectomy that directly impact patient outcomes

Navigation & Access

Difficult or tortuous anatomy, small vessel size, or heavily calcified arteries can make access difficult

Incomplete Clot Removal

Some devices may not capture or aspirate the entire thrombus, especially if the clot is old, organized, or adherent to the vessel wall

Distal Embolization

During clot retrieval, parts of the thrombus can break off and travel downstream to distal or previously unaffected territories

Risk of Vessel Injury

Manipulation in delicate or diseased vessels can cause dissection or perforation

Time is Brain: Pressure

Every minute without reperfusion leads to neuron loss

FIRST PASS EFFECT (FPE) → BETTER OUTCOMES

FPE (full recanalization in the 1st Pass) is the most powerful predictor of clinical outcome with best safety results

NASA registry (25 U.S. centers)	First Pass Effect Group (n=89)	Total Patient Population (n=354)	→ FPE: 25.1% of total patient population
MCA occlusions	64.0%	52.5%	→ FPE: more commonly observed in MCA occlusions
ICA occlusions	10.1%	27.7%	
BGCs use	64%	34.7%	→ FPE: more commonly observed when balloon guide catheters were used
Median time to revascularization	34 min	60 min	→ Median time to revascularization was significantly faster in FPE group (p=.0003)
mRS ≤ 2	61.3%	35.3%	→ Patient outcomes were significantly better in FPE group (p=.013)



FIRST PASS EFFECT (FPE) → BETTER OUTCOMES

FPE (full recanalization in the 1st Pass) is the most powerful predictor of clinical outcome with best safety results

Faster Recanalization

“Time is brain”


Lower Complication Rates

Thanks to smaller number of passes


Complete Reperfusion

Better quality recanalization results in better outcomes


Each hour of delay to treatment is associated with a 5.5% absolute decline in the likelihood of good outcome



Proportion of good outcome by mTICI grade (p<0.0001 for overall comparison)



CAUSES OF RECANALIZATION FAILURE




An organized (hard, fibrin-rich) clot is more resilient and less sticky than a fresh (soft, red blood cell-rich) clot.

As a result, it engages less effectively with a typical stent retriever, increasing the risk of clot disengagement or “missed” retrieval—particularly in the setting of a tortuous arterial anatomy.

CAUSES OF RECANALIZATION FAILURE

Up to 40% of thrombectomies result in clot particle escape to distal or previously unaffected vascular territories, adversely affecting outcomes



Aspiration as well as open-ended stent retrievers do not protect against the risk of clot escape to distal vasculature

DESIGN

Designed for 1st PASS SUCCESS with ALL Clot Types



NevaTM

NOT JUST ANOTHER STENT-RETRIEVER



All-Clots Capability

Remove even the hardest, fibrin-rich clots

Contain soft, fragment-prone clots that easily disintegrate

First-Pass Success

Reduce time-to-recanalization

Improve on recanalization quality, with higher TICI 2c/3 rates


Ease of Use

Allow real-time feedback during retrieval

Synergistic with all access philosophies


TARGETING FIRST PASS SUCCESS WITH ALL CLOT TYPES

NevaTM DESIGNED WITH DROP ZONETM TECHNOLOGY TO
CAPTURE ALL TYPES OF CLOT INSIDE



Drop Zones offset
at 90°, acting as
pockets to capture
thrombi inside

NevaTM **NET** THE FIRST & ONLY MICRO-FILTRATION TECHNOLOGY DESIGNED TO
MAXIMIZE CLOT RETENTION



32 intricately braided
nitinol strands of
.00125" creating a filter
with an average pore
size of 385.3±68 µm

TARGETING FIRST PASS SUCCESS WITH ALL CLOT TYPES

DROP ZONES™

Offset at 90°, act as entry points to laterally integrate all clot types for fast and effective recanalization

CAREFULLY CALIBRATED RADIAL FORCE

for wall apposition and clot retention during retrieval

LARGE OPENINGS AND CLOSED ENDS

for endothelial tissue protection

CLOSED DISTAL TIP

Clot gets inside,
Clot stays inside!


RADIOPAQUE MARKERS

for positioning & clot engagement guidance


Inside the NeVa NET distal tip, 32 intricately braided strands create a **385.3±68 µm PORE-SIZED FILTER**

NEVA DESIGN EVOLVED WITH FEEDBACK WITH AN AMBITION TO EXCEED AND EXCEL

2018




2024



- Flow restoration zone exchanged for an additional Drop Zone
- Several sizes piloted, including an ultra-long 5 Drop Zone NeVa
- Addition of “zebra” fluoroscopy markers to the pusher-wire
- Lower profile markers (maintained visibility) ... and many minor design adjustments to improve microcatheter compatibility, navigability and other features
- NET: the first ever integrated microfilter in the basket structure

NEVA NET HAS A SINGULAR MICRO FILTRATION TECHNOLOGY FOR MAXIMAL CLOT RETENTION




**Neva[™] NET**

32 intricately
braided nitinol
strands of .00125"
creating a filter with
an average pore size
of $385.3 \pm 68 \mu\text{m}$

NEVA IS DESIGNED FOR WALL APPPOSITION

Compressive Radial Force Measurements



NEVA 3 MM: DROP ZONE TECHNOLOGY REACHING NEW ANATOMY

**A SEAMLESS,
CONTINUOUS STRUCTURE**
gliding with ease in
.017" ID microcatheters


**FOR ARTERIES
AS SMALL AS 1.5 MM:**
perfect for
M2, ACA & PCA occlusions

Radial force **ATTUNED TO
DELICATE ANATOMY**



Extending the reach of Vesalio's proven Drop Zone™ technology engineered for **ALL CLOT TYPES**

NEVA 3MM: DROP ZONE TECHNOLOGY REACHING NEW ANATOMY



Pusherwire Specifications

- Diameter: 0.014", .017" Microcatheter compatible
- Similar tapering, miniaturized
- Zebra markers: 37.5 cm from the proximal end of the pusherwire

Basket Specifications

- Reduced radial force
- Smooth, continuous structure
- 2 Drop Zones
- Indicated in 1.5 - 3 mm arteries

CLINICAL DATA & EXPERIENCE

NevaTM

Designed for 1st PASS SUCCESS with ALL Clot Types






LV-MKT-010 REV E



CLINICAL EXPERIENCE

- More Than 15000 Cases Worldwide
- 16+ publications
- All-Clots Capability
- Superior First-Pass Success
- Maximized Clot Retention
- Proven Safety



Neva Data

21

Excalibur

FPE 94%

Wield Your Weapon Wisely!!!

Copyright © Brain. All rights reserved.



NEVA PRE-CLINICAL TESTING



PRE-LAUNCH ANIMAL TESTING: ALL CLOTS CAPABILITY PROVEN

Clot Type	Soft	Hard	Ultra Hard	All Clots
Clot morphology	Whole Blood "RED" Clot	Plasma Rich "WHITE" Clot	Clot modeled from ONYX 500	RED, WHITE and ONYX 500
N =	19	5	11	35
Length of clots - mm	10-40	6-12	4-12	4-40
1 st Pass TIC1 3	84%	60%	55%	71%
Final TIC1 3	89%	NR	82%	83%
Final TIC1 2b/3	100%	100%	91%	97%
Average # of passes for final recanalization	1,05	1,00	1,63	1,23



REMOVED
ALL
SOFT AND
HARD CLOTS IN
A SINGLE PASS

REMOVED
91%
OF THE ULTRA
HARD CLOTS

ACHIEVED
71% 97%
FIRST-PASS FINAL
RECANALIZATION ACROSS
ALL CLOT TYPES





PRE-LAUNCH IN VITRO TESTING:
CONSISTENT EFFECTIVENESS AT REMOVING ORGANIZED CLOTS

Benchtop testing comparing all available devices in 2016¹ ***“All stent retrievers failed when interacting with large white thrombi (≥ 6mm)”***

Same study design, same model, same physician
DIFFERENT RESULTS in 2018²

Solitaire*:	0/5	Trevo:	0/5
Embotrap*:	0/5	Eric:	0/5
Preset*:	0/5	Preset LT:	0/5
Catch*:	0/5	Separator 3D:	0/5
Revive*:	0/5	Mindframe:	0/5

NeVa: 6/10 successful complete removals of large (≥ 6 mm), white (hard) thrombi



NEVA CLINICAL EXPERIENCE






CLINICAL EXPERIENCE IN FIRST-LINE SUPERIOR FIRST-PASS SUCCESS

First-Pass Recanalization

■ TICI 2C-3 ■ TICI 2B-3




1. Zaidat OO et al., Primary Results of the Multicenter ARISE II Study (Analysis of Revascularization in Ischemic Stroke With EmboTrap), *Stroke* Volume 49, Number 5
2. Nogueira RG et al., Thrombectomy With the pRESET vs Solitaire Stent Retrievers as First-Line Large Vessel Occlusion Stroke Treatment: A Randomized Clinical Trial. *JAMA Neurol.* 2024 Feb 1;81(2):170-178.
3. Abbasi et al., *JNS* 2021;13:212-216. 2. Yoo AJ, et al., *JNS* 2023. doi: 10.1136/jnis-2023-020960.
4. Yoo AJ, Geyik S, Froehler MT, et al., Primary results from the CLEAR study of a novel stent retriever with drop zone technology, *Journal of NeuroInterventional Surgery* 2024;16:1220-1227
5. Bajrami A, Ertugrul O, Senadim S, et al. First-pass results of mechanical thrombectomy with the two-drop zone NeVa™ device. *Interv Neuroradiol.* 2022;28:255-263.





CLINICAL EXPERIENCE IN FIRST-LINE RECANALIZATION BEFORE RESCUE

Recanalization after 2-3 passes – without rescue




1. Zaidat OO et al., Primary Results of the Multicenter ARISE II Study (Analysis of Revascularization in Ischemic Stroke With EmboTrap), *Stroke* Volume 49, Number 5
2. Nogueira RG et al., Thrombectomy With the pRESET vs Solitaire Stent Retrievers as First-Line Large Vessel Occlusion Stroke Treatment: A Randomized Clinical Trial. *JAMA Neurol.* 2024 Feb 1;81(2):170-178.
3. Yoo AJ, Geyik S, Froehler MT, et al., Primary results from the CLEAR study of a novel stent retriever with drop zone technology, *Journal of NeuroInterventional Surgery* 2024;16:1220-1227
4. Bajrami A, Ertugrul O, Senadim S, et al. First-pass results of mechanical thrombectomy with the two-drop zone NeVa™ device. *Interv Neuroradiol.* 2022;28:255-263.





CLINICAL EXPERIENCE IN FIRST-LINE FINAL RECANALIZATION

Final Recanalization Rates



1. Zaidat OO et al., Primary Results of the Multicenter ARISE II Study (Analysis of Revascularization in Ischemic Stroke With EmboTrap), *Stroke* Volume 49, Number 5
2. Nogueira RG et al., Thrombectomy With the pRESET vs Solitaire Stent Retrievers as First-Line Large Vessel Occlusion Stroke Treatment: A Randomized Clinical Trial. *JAMA Neurol.* 2024 Feb 1;81(2):170-178.
3. Yoo AJ, Geyik S, Froehler MT, et al., Primary results from the CLEAR study of a novel stent retriever with drop zone technology, *Journal of NeuroInterventional Surgery* 2024;16:1220-1227
4. Bajrami A, Ertugrul O, Senadim S, et al. First-pass results of mechanical thrombectomy with the two-drop zone NeVa™ device. *Interv Neuroradiol.* 2022;28:255-263.






CLINICAL EXPERIENCE IN RESCUE SIGNIFICANT IMPACT

Significant improvement was observed with NeVa as a rescue device after failed initial attempts

Recanalization Improvement



CLEAR STUDY RESULTS IN DETAIL







CLEAR STUDY EFFICACY ENDPOINT (mITT, n=107) FIRST-PASS REPERFUSION

NeVa achieved a superior rate of first pass successful reperfusion compared to predicate devices - prespecified secondary endpoint hierarchical testing

SUCCESSFUL RECANALIZATION % eTICI 2B-3



EXCELLENT RECANALIZATION % eTICI 2C-3





CLEAR STUDY EFFICACY ENDPOINT (mITT, n=107) FIRST-PASS REPERFUSION

NeVa achieved a superior rate of first pass successful reperfusion compared to predicate devices - prespecified secondary endpoint hierarchical testing

SUCCESSFUL RECANALIZATION
% eTICI 2B-3

≈ 7.5
/
10


EXCELLENT RECANALIZATION
% eTICI 2C-3

≈ 5
/
10





CLEAR STUDY EFFICACY ENDPOINT (mITT, n=107) eTICI 2B-3 WITHIN 3 PASSES




NeVa successfully demonstrated non-inferiority to the performance goal based on predicate devices (prespecified primary efficacy analysis)

NeVa further demonstrated post hoc superiority to the predicate performance goal





CLEAR STUDY SAFETY ENDPOINT (ITT, n=139) COMPOSITE 90-DAY MORTALITY & 24-HR sICH



Composite 90-DAY all-cause mortality and/or 24-hour sICH met the performance goal (non-significantly lower)






CLEAR STUDY CLINICAL OUTCOMES

90-DAY mRS & OTHER

90-DAY MRS ≤ 2 (mITT, n=106)

Outcomes with NeVa were superior to predicate studies in exploratory analysis

65.1% vs. PG: 54.7%, p = 0.03



OTHER CLINICAL OUTCOMES (mITT, n=107)

NIHSS at 24-hours, median (IQR) 4 (1-8)

NIHSS change from baseline to 24-hrs, median (IQR) -10 (-15 to -5)

NIHSS at 5-10 days or discharge, median (IQR) 2 (0-5)


NIHSS change from baseline to 5-10 days or discharge, median (IQR) -11 (-17 to -8)




NEVA NET CLINICAL EXPERIENCE




NUMBER OF FRAGMENTS
(0.2 – 1 MM)¹




NUMBER OF
FRAGMENTS (> 1 MM)¹



Median number of
clot fragments
generated by NeVa
NET were
1/4th
those generated by
Solitaire¹



RETENTION OF FRAGMENT SIZES > 1 MM²



NeVa NET was
significantly better
at preventing clot
fragments >1mm
from embolizing
distal territories²

The injury scores were nearly identical to predicate and primarily related to endothelial loss, occasional IEL disruption and limited medial injury

- NeVa NET 5.5mm compared to Solitaire 4x40 in vessels < 3mm
Solitaire 6x40 in vessels > 3mm
- Soft and hard platelet rich thrombus were delivered in 6 vascular territories
- 4 retrievals performed in each vessel
- Vessel sizes ranged from 1.8 to 3.4 mm
- Harvested vessels analyzed by an independent veterinary pathologist and compared for thrombectomy induced acute vascular injury

Example

B. SolitaireTM 4x40

MEDIA COMPRESSIVE HYPOCELLULARITY/NECROSIS SCORE 2 | ENDOTHELIAL LOSS SCORE 4




Figure 2. Animal 7173 (H&E). Left Internal Carotid. Solitaire 4x40mm (hard clot). Day 0 (acute). The media demonstrating pale eosinophilic hypocellularity and scattered karyorrhexis (black arrowheads) characteristic of media necrosis and cell death; M, area bound by black dotted lines = area of normal media; clear arrowheads = complete endothelial cell loss/denudation; A = adventitia.

A. NeVa NETTM 5.5x30

MEDIA COMPRESSIVE HYPOCELLULARITY/NECROSIS SCORE 2 | ENDOTHELIAL LOSS SCORE 4




Figure 3. Animal 7173 (H&E). Right Internal Carotid. NeVa5.5 Net (hard clot). Day 0 (acute). Partial thickness pale eosinophilic hypocellularity (asterisks) and scattered karyorrhexis (black arrowheads) within the media (M); clear arrowheads = complete endothelial cell loss/denudation; A = adventitia.

- A = adventitia,
- M = media
- ▶ = scattered hypocellularity and karyorrhexis
- ▽ = endothelial loss
- * = pale eosinophilic hypocellularity



Schwab et al. JNIS

1st publication on NeVa NET 5.5 x 37 mm

- 51 patients
- 3 EU sites
- >70% of patients had ICA occlusions with extensive clot burden (mean: 25.1±13.3 mm, range 4-50 mm)

Bhogal et al. Acta Neurochirurgica

1st Publication on NeVa NET 4.0 x 30 mm

- 20 patients
- 2 EU sites
- Mixed occlusion locations with 50% M1 and 40% ICA

1ST PASS RATES
HIGH

eTICI 2c-3
55 %

eTICI 2b-3
78.5 %

eTICI 2c-3
94 %

EMBOLIZATION TO
NEW TERRITORIES
LOW

3.9 %

0

REFERENCES

1. Ulm A, Khachatryan T, Grigorian A, Nogueira R. Preclinical evaluation of the NeVa™ stent retriever: safety and efficacy in the swine thrombectomy model. *Intervent Neurol.* 2018;7:205-217.
2. Machi P, Ulm A, Bernava G, et al. Experimental evaluation of the NeVa™ thrombectomy device: a novel stent retriever conceived to improve the efficacy of organized clot removal. *J Neuroradiol.* 2019;46:163-167.
3. Ribo M, Requena M, Macho J, et al. Mechanical thrombectomy with a novel stent retriever with multi-functional zones: initial clinical experience with the NeVa™ thrombectomy device. *J Neuroradiol.* 2020;47:301-305.
4. Akpınar CK, Ozdemir AO, Gurkas E, et al. Favorable first-pass recanalization rates with the NeVa™ thrombectomy device in acute stroke patients: initial clinical experience. *Interv Neuroradiol.* 2020;26:1-7.
5. Borggrefe J, Goertz L, Abdullayev N, et al. Mechanical thrombectomy with the novel NeVa™ M1 stent retriever: Do the drop zones represent a risk or benefit? *World Neurosurgery.* 2021;150:e334-e342.
6. Bajrami A, Ertugrul O, Senadim S, et al. First-pass results of mechanical thrombectomy with the two-drop zone NeVa™ device. *Interv Neuroradiol.* 2022;28:255-263.
7. Masthoff M, Krähling H, Akkurt BH, et al. Evaluation of the effectiveness and safety of the multizone NeVa™ stent retriever for mechanical thrombectomy in ischemic stroke. *Neuroradiology.* 2023;65:1-10
8. Melki G-J, Demeestere J, Laenen A et al., Immediate and 90-Day Clinical Outcome of Patients with Acute Stroke Treated with the NeVa-Vesalio Mechanical Thrombectomy Device: A Retrospective Case Series, *World Neurosurgery*, 2023, Pages e212-e221
9. Habibi, MA et al., The safety and efficacy of NeVa mechanical thrombectomy device in acute ischemic stroke: A systematic review and meta-analysis, *Journal of Clinical Neuroscience*, Volume 130, 110892
10. Yoo AJ, Geyik S, Froehler MT, et al., Primary results from the CLEAR study of a novel stent retriever with drop zone technology, *Journal of NeuroInterventional Surgery* 2024;16:1220-1227
11. Bhogal P, Mancuso-Marcello M, Fairhead R. et al. The NeVa stent-retriever - a single-centre real-world experience. *Interv Neuroradiol.* 2025 May 21:15910199251337176. doi: 10.1177/15910199251337176. Epub ahead of print. PMID: 40398472; PMCID: PMC12095227
12. Fairhead R, Mancuso-Marcello M, Ahmed R, et al., The use of the NeVa stent-retriever for bail-out mechanical thrombectomy. *Interv Neuroradiol.* 2025 Sep 29:15910199251381491. doi: 10.1177/15910199251381491. Epub ahead of print. PMID: 41021792; PMCID: PMC12479453
13. Li et al., Impact of stent-retriever tip design on distal embolization during mechanical thrombectomy: a randomized in vitro evaluation, *JNIS*, May 2023
14. Anagnostakou et al., Preclinical safety and efficacy of the NeVa NET™: A novel thrombectomy device with integrated embolic distal protection, *JVIN*, Feb 2023
15. Schwab R, Kis B, Réka BA, et al. First clinical multicenter experience of the new NeVa NET 5.5 thrombectomy device, *J NeuroIntervent Surg* Epub ahead of print: 16 June 2025. doi:10.1136/jnis-2025023476
16. Bhogal P, et al., The NeVa Net stent-retriever – initial report of 20 cases from two high volume centres, *Acta Neurochirurgica*, July 2025



DO THE DROP ZONE

Recommendations for new NeVa and NeVa NET users

NeVaTM

Designed for 1st PASS SUCCESS with ALL Clot Types




NEVA DESIGN, TIPS & TRICKS

- **Features of NeVa**
- **Choosing The Correct Neva Size**
- **Positioning Neva**
- **Deploying Neva**
- **Retrieving Neva**

DROP ZONE™ THE CLOT INSIDE




DESIGNED WITH DROP ZONE™ TECHNOLOGY TO
CAPTURE ALL TYPES OF CLOT INSIDE



Drop Zones offset
at 90°, acting as
pockets to capture
thrombi inside




THE FIRST & ONLY MICRO-FILTRATION TECHNOLOGY DESIGNED TO
MAXIMIZE CLOT RETENTION



32 intricately braided
nitinol strands of
.00125" creating a filter
with an average pore
size of 385.3±68 µm



Neva™ & Neva™ NET STRUCTURE



NeVa™ & NeVa™ NET SIZES

	Product Number (SKU)	Maximal Diameter	Working Length (mm)	# of Drop Zones	Total Basket Length (mm)	Min MC ID	Pusher wire Length (cm)	Recommended Vessel Diameter
NeVa 3.0 x 17	VN-3017-32RR	30	17	2	30	0.017	200	≥ 1.5 and ≤ 3.0
NeVa 4.0 x 22	30020V-MS	40	22	2	39	0.021	200	≥ 2.0 and ≤ 3.5
NeVa 4.5 x 29	VN-4529-03RR	45	29	3	46	0.021	200	≥ 2.0 and ≤ 4.5
NeVa 5.5 x 37	VN-5537-03RR	55	37	3	56	0.027	200	≥ 3.5 and ≤ 5.5
NeVa NET 4.0 x 30	VN-4030-03NR	40	30	3	49	0.021	200	≥ 2.0 and ≤ 3.5
NeVa NET 5.5 x 37	VN-5537-03NR	55	37	3	56	0.027	200	≥ 3.5 and ≤ 5.5






NEVA DESIGN, TIPS & TRICKS

- Features of NeVa
- **Choosing The Correct Neva Size**
- Positioning Neva
- Deploying Neva
- Retrieving Neva


NeVa™ 5.5

Ideal for proximal occlusions all the way into proximal M1



Neva™ 3.0

Use in the MCA (M2), the ACA, as well as similar caliber arteries in the posterior circulation, such as the PCA




NEVA DESIGN, TIPS & TRICKS

- Features of NeVa
- Choosing The Correct Neva Size
- **Positioning Neva**
- Deploying Neva
- Retrieving Neva



Position distal enough to allow all Drop Zones to interact with clot


Balance the benefit & risk of distal placement




Position distal enough to allow all Drop Zones to interact with clot

Balance the benefit & risk of distal placement

Place the proximal marker at the edge of the occlusion





Place the proximal marker at the edge of the occlusion



All Drop Zones are interacting with the clot




Positioning not distal enough



The first Drop Zone is deployed proximal to clot





To place the proximal marker at the edge of the occlusion you need to go one full NeVa length beyond the occlusion site

	Drop Zones	Full Length	
NeVa 3.0 x 17	2	30 mm	≈ 3 cm
NeVa 4.0 x 22	2	39 mm	≈ 4 cm
NeVa NET 4.0 x 30	3	49 mm	≈ 5 cm
NeVa 4.5 x 29	3	46 mm	≈ 4.5 cm
NeVa 5.5 x 37	3	56 mm	≈ 5.5 cm
NeVa NET 5.5 x 37	3	56 mm	≈ 5.5 cm



MICROCATHETER CONSIDERATIONS

	Drop Zones	Full Length	MC Compatibility	
NeVa 3.0 x 17	2	30 mm	0.017"	<ol style="list-style-type: none">1. Always flush NeVa before insertion2. Choose micro-catheters with good distal support3. In tortuous anatomies, consider using larger microcatheters for better distal support
NeVa 4.0 x 22	2	39 mm		
NeVa NET 4.0 x 30	3	49 mm	0.021"	
NeVa 4.5 x 29	3	46 mm		
NeVa 5.5 x 37	3	56 mm	0.027"	
NeVa NET 5.5 x 37	3	56 mm		



NEVA DESIGN, TIPS & TRICKS

- Features of NeVa
- Choosing The Correct Neva Size
- Positioning Neva
- **Deploying Neva**
- Retrieving Neva



Release tension on the micro-catheter just before starting the deployment

Deploy slowly and carefully to protect your ideal position

Expect initial anchoring after 1 cm of unsheathing



NEVA DESIGN, TIPS & TRICKS

- Features of NeVa
- Choosing The Correct Neva Size
- Positioning Neva
- Deploying Neva
- **Retrieving Neva**

NeVa™ RETRIEVAL

SHOULD YOU WAIT ?



You can but do not have to wait


Check the Drop Zone markers:

Are they fully open?

Is the proximal marker at the edge of the occlusion?

Check the position of your DAC if you're using one

START SLOW PULL



Apply slow & gentle vessel straightening traction to start pulling NeVa proximally

Start co-aspiration from your proximal catheter

WATCH THE DROP ZONE MARKERS




Watch the Drop Zone markers, observe if one of the pairs is compressing on one another


Neva™ RETRIEVAL

MARKERS COMPRESSED TOGETHER

you may be
adjacent to a hard clot




SLOW DOWN!




MARKERS SPRING OPEN

you may now be
at the proximal edge of the hard clot




THE DROP ZONE IS ON THE CLOT



NeVa™ ANGLE OF RETRIEVAL




OPTION 1:
Bring NeVa proximally towards the DAC and align NeVa with the tip of your DAC



OPTION 2:
Use NeVa as an anchor, and drive up your DAC

1. Straightening the anatomy eases retrieval
2. Avoids clot fragmentation
3. Aspiration via DAC will be more efficient






RETRIEVING LARGE, HARD CLOTS

After deploying NeVa, bring the DAC tip up to the proximal marker



Remove excess tension from the DAC and slowly retrieve NeVa.

If significant resistance is encountered, stop retrieval to avoid shearing. Clot may be partially incorporated and trapped between stent and DAC.



Tighten the RHV of DAC around the MC and retrieve the whole system together (DAC + MC + NeVa) while gently aspirating

