

Neva™



DESIGNED FOR  
**FIRST-PASS SUCCESS**  
WITH **ALL CLOT TYPES**




Neva  **NET** DESIGNED TO  
**MAXIMIZE CLOT RETENTION**



# Neva NET™

DESIGNED TO

**MAXIMIZE CLOT RETENTION**



**“A first pass rate in a league of its own...”**

Dr Paul Bhogal, UK

**Exceptional first-pass eTICI 2c-3 rates:**

**94%**

in a standard stroke  
practice cohort<sup>1</sup>

**55%**

in a challenging proximal  
dominated cohort (>70% ICA)<sup>2</sup>

**“A lower distal embolization rate,  
especially in hard clots...”**

Dr Roland Schwab, Germany

**Reduced fragment loss versus competitive stent retrievers:**

**78%** reduction in small fragments  
(0.2 – 1 mm) vs SR1\* (p = .048)<sup>3</sup>

\*SR1: Solitaire Ref: 3

# Neva™ 3mm

DROP ZONE™ TECHNOLOGY


## REACHING NEW ANATOMY

Extending the reach of Vesalio's proven Drop Zone™ technology engineered for All Clot Types

For arteries as small as 1.5 mm  
**Perfect for M2, ACA & PCA occlusions**

Radial force  
**attuned to delicate anatomy**

A seamless, continuous structure  
**gliding with ease in .017" ID microcatheters**





# SUPERIOR FIRST-PASS SUCCESS, PROVEN SAFETY<sup>1</sup>

DROP ZONE TECHNOLOGY™ CAPTURE CLOT INSIDE

Offset at 90° **Drop Zones** act as entry points to laterally integrate all clot types for fast and effective recanalization

**Closed Distal Tip**  
Clot gets inside  
Clot stays inside

**Large openings and closed ends** for endothelial tissue protection

**Carefully calibrated radial force** for wall apposition and clot retention during retrieval

Inside the **NeVa NET** distal tip, 32 intricately braided nitinol strands create a **385.3±68 µm pore sized filter**

**Radiopaque Markers** for positioning & clot engagement guidance

Drop Zones remove

**ALL CLOT TYPES**

including ultra-hard, long ones<sup>5-6</sup>

NeVa NET first-pass  
TICI 2c-3 up to:

**94%**

in early studies<sup>1</sup>

NeVa first-pass TICI  
2c-3 up to:

**67%**

in published studies  
Range: 45% - 67%<sup>7-14</sup>

Recanalization within  
3 NeVa passes, with no  
rescue, up to:

**92%**

in published studies  
Range: 60% - 92%<sup>7-14</sup>

# REFERENCES

1. Bhogal P, et al. The NeVa Net stent-retriever – initial report of 20 cases from two high-volume centres. *Acta Neurochirurgica*. 2025 Jul.
2. Schwab R, et al. First clinical multicenter experience of the new NeVa NET 5.5 thrombectomy device. *JNIS*. 2025 Jun 16. [Epub ahead of print]
3. Anagnostakou, et al. Preclinical safety and efficacy of the NeVa NET™: A novel thrombectomy device with integrated embolic distal protection. *JVIN*. 2023 Feb.
4. Li, et al. Impact of stent-retriever tip design on distal embolization during mechanical thrombectomy: a randomized in vitro evaluation. *JNIS*. 2023 May.
5. Ulm A, et al. Preclinical evaluation of the NeVa™ stent retriever: safety and efficacy in the swine thrombectomy model. *Intervent Neurol*. 2018;7:205–217.
6. Machi P, et al. Experimental evaluation of the NeVa™ thrombectomy device: a novel stent retriever conceived to improve the efficacy of organized clot removal. *J Neuroradiol*. 2019;46:163–167.
7. Ribo M, et al. Mechanical thrombectomy with a novel stent retriever with multi-functional zones: initial clinical experience with the NeVa™ thrombectomy device. *J Neuroradiol*. 2020;47:301–305.
8. Akpınar CK, et al. Favorable first-pass recanalization rates with the NeVa™ thrombectomy device in acute stroke patients: initial clinical experience. *Interv Neurol*. 2020;26:1–7.
9. Borggrefe J, et al. Mechanical thrombectomy with the novel NeVa™ M1 stent retriever: Do the drop zones represent a risk or benefit? *World Neurosurgery*. 2021;150:e334–e342.
10. Bajrami A, et al. First-pass results of mechanical thrombectomy with the two-drop zone NeVa™ device. *Interv Neurol*. 2022;28:255–263.
11. Masthoff M, et al. Evaluation of the effectiveness and safety of the multizone NeVa™ stent retriever for mechanical thrombectomy in ischemic stroke. *Neuroradiology*. 2023;65:1–10.
12. Melki G-J, et al. Immediate and 90-day clinical outcome of patients with acute stroke treated with the NeVa-Vesalio mechanical thrombectomy device: A retrospective case series. *World Neurosurgery*. 2023; e212–e221.
13. Yoo AJ, et al. Primary results from the CLEAR study of a novel stent retriever with drop zone technology. *J NeuroInterv Surg*. 2024;16:1220–1227.
14. Bhogal P, et al. The NeVa stent-retriever – a single-centre real-world experience. *Interv Neurol*. 2025 May 21;15910199251337176. doi: 10.1177/15910199251337176. Epub ahead of print. PMID: 40398472; PMCID: PMC12095227.

Prescription use only. Please see [www.vesalio.com](http://www.vesalio.com) for indications for use, contraindications, warnings, precautions, and potential adverse events. For complete information, including detailed instructions for use, please refer to the Instructions for use.


©Vesalio, Inc. 2026. All rights reserved. The name Vesalio and the product names and associated logos are registered and unregistered trademarks of Vesalio, Inc. The following trademarks are registered at the USPTO: Vesalio®, NeVa®, NeVa VS®, NeVa NET®, and pVasc®. The following trademarks are registered at the EUIPO and UKIPO: Vesalio™, Drop Zone™, NeVa™, pVasc™, and enVast™.

Product Name	Product Code	Maximal Diameter (mm)	Working Length (mm)	Full Length (mm)	Nr. of Drop Zones	Recommended Vessel diameter (3mm)	Min. MC ID	Pusher Wire (cm)
NeVa 3.0 x 17	VN-3017-32RR	3.0	17	30	2	≥ 1.5 and ≤ 3.0	0.017"	200
NeVa 4.0 x 22	30020V-MS	4.0	22	39	2	≥ 2.0 and ≤ 3.5	0.021"	200
NeVa 4.5 x 29	VN-4529-03RR	4.5	29	46	3	≥ 2.0 and ≤ 4.5	0.021"	200
NeVa 5.5 x 37	VN-5537-03RR	5.5	37	56	3	≥ 3.5 and ≤ 5.5	0.027"	200
NeVa NET 4.0 x 30	VN-4030-03NR	4.0	30	49	3	≥ 2.0 and ≤ 3.5	0.021"	200
NeVa NET 5.5 x 37	VN-5537-03NR	5.5	37	56	3	≥ 3.5 and ≤ 5.5	0.027"	200

# SUPERIOR **FIRST PASS** SUCCESS

## SUCCESSFUL FIRST-PASS RECANALIZATION eTICI 2B-3

(mITT, n = 107, Superiority on the Performance Goal, p = 0.0008)




## SUCCESSFUL FIRST-PASS RECANALIZATION eTICI 2C-3

(mITT, n = 107, Numerically higher than both predicates: ARISE 2: 40%, TIGER: 41.4%)




## COMPOSITE PRIMARY SAFETY ANALYSIS

90-DAY ALL-CAUSE MORTALITY AND /OR POST PROCEDURE sICH (ITT, n = 139)




## CLINICAL OUTCOMES

90-DAY MRS 0-2 (mITT, n = 106, vs. performance goal: 54.7%, p = 0.03)



Yoo AJ, et al. Primary results from the CLEAR study of a novel stent retriever with drop zone technology, JNIS, 2023



# MAXIMIZE CLOT RETENTION

**SIGNIFICANTLY BETTER AT PREVENTING CLOT FRAGMENTS FROM EMBOLIZING DISTAL TERRITORIES<sup>1,2</sup>**

**28% – 63% LESS**

FRAGMENTS > 1 MM  
GENERATED VS. SOLITAIRE  
P = .003<sup>1</sup> ; P = .037<sup>2</sup>

**24% LESS**

FRAGMENTS > 1 MM  
GENERATED VS. EMBOTRAP  
P = .003<sup>1</sup>

**78% LESS**

FRAGMENTS 0.2 - 1 MM  
GENERATED VS. SOLITAIRE  
P = .0048<sup>2</sup>

## EXCEPTIONAL RECANALIZATION SUCCESS<sup>3,4</sup>

### RECANALIZATION RATES FROM INITIAL MULTI-CENTER EXPERIENCE

#### NEVA NET 5.5 MM

51 PATIENTS FROM 2 EUROPEAN SITES<sup>3</sup>:

- 71% ICA occlusions
- Mean Clot Length: 25.1 +/- 13.3 mm (range 4–50 mm)

**55%**

FIRST-PASS  
TICI 2C-3

**79%**

FIRST-PASS  
TICI 2B-3

**98%**

FINAL  
TICI 2B-3

#### NEVA NET 4.0 MM

20 PATIENTS FROM 2 EUROPEAN SITES<sup>4</sup>:

- 50% M1, 40% ICA occlusions
- 70% with hyperdense clot
- Mean Clot Length: 21.6 +/- 8.9 mm

**94%**

FIRST-PASS (& FINAL)  
TICI 2C-3

1.Li J, et al. Impact of stent-retriever tip design on distal embolization during mechanical thrombectomy: a randomized in vitro evaluation, JNIS, 2023


2.Anagnostakou V, et al. Preclinical safety and efficacy of the NeVa NET™: A novel thrombectomy device with integrated embolic distal protection: Preclinical safety and efficacy of the NeVa NET™. JVIN, 2023

3.Schwab et al., First clinical multicenter experience of the new NeVa NET 5.5 thrombectomy device, J NeuroIntervent Surg, 2025

4.Bhogal P. et al., The NeVa Net stent-retriever – initial report of 20 cases from two high volume centres, Acta Neurochirurgica, 2025

## EFFECTIVE IN RETRIEVING ALL CLOT TYPES

### NEVA SUCCEEDS WITH ALL CLOT TYPES<sup>1</sup>




### NEVA SUCCEEDS WHERE OTHER STENT RETRIEVERS FAIL<sup>2,3</sup>

ALL TESTED DEVICES FAILED IN REMOVING LARGE (> 6 MM), WHITE THROMBI<sup>2</sup>

NEVA SUCCEEDED WHERE ALL OTHERS FAILED<sup>3</sup>

Solitaire\*  
Trepo  
Embotrap\*  
Preset\*  
Preset LT  
Catch\*  
Eric  
Separator 3D  
...\*\*

0%



60%

\* These devices were able to minimally displace clots without removing them,

\*\* The list of devices that failed is incomplete

1.Ulm A, et al. Preclinical evaluation of the NeVa™ stent retriever: safety and efficacy in the swine thrombectomy model, Intervent Neurol, 2018


2.Machi P, et al. Experimental evaluation of stent retrievers' mechanical properties and effectiveness, JNIS, 2016

3.Machi P, et al. Experimental evaluation of the NeVa™ thrombectomy device a novel stent retriever conceived to improve efficacy of organized clot removal. J Neuroradiology, 2019




# EFFECTIVE RECANALIZATION IN FIRST LINE USE


## FIRST-PASS RECANALIZATION RATES



## RECANALIZATION AFTER 2-3 PASSES – WITHOUT RESCUE



## FINAL RECANALIZATION RATES





1.Zaidat OO et al., Primary Results of the Multicenter ARISE II Study (Analysis of Revascularization in Ischemic Stroke With EmboTrap), Stroke Volume 49, Number 5  
 2.Nogueira RG et al., Thrombectomy With the pRESET vs Solitaire Stent Retrievers as First-Line Large Vessel Occlusion Stroke Treatment: A Randomized Clinical Trial. JAMA Neurol. 2024 Feb 1;81(2):170-178.  
 3.Yoo AJ, Geyik S, Froehler MT, et al., Primary results from the CLEAR study of a novel stent retriever with drop zone technology. Journal of NeuroInterventional Surgery 2024;16:1220-1227  
 4.Bajrami A, Ertugrul O, Senadim S, et al. First-pass results of mechanical thrombectomy with the two-drop zone NeVa™ device. Interv Neuroradiol. 2022;28:255-263  
 5.Abbasi et al., JNS 2021;13:212-216, 2. Yoo AJ, et al., JNS 2023. doi:10.1136/jnis-2023-020960.



# SIGNIFICANT IMPACT IN RESCUE

SIGNIFICANT IMPROVEMENT WAS OBSERVED WITH NEVA AS A RESCUE DEVICE AFTER FAILED INITIAL ATTEMPTS



<sup>1</sup>Masthoff M, et al., Evaluation of the effectiveness and safety of the multizone NeVa™ stent retriever for mechanical thrombectomy in ischemic stroke. *Neuroradiology*. 2023;65:1-10  
<sup>2</sup>Fairhead R, et al., The use of the NeVa stent-retriever for bail-out mechanical thrombectomy. *Interv Neuroradiol*. 2025 Sep