


DROP ZONE™ the CLOT INSIDE

DROP ZONE TECHNOLOGY

taking clot capture to the next level



BALANCED DESIGN
optimized radial force
balanced with large openings

SMART MARKERS
real-time feedback
during retrieval

CLOSED DISTAL TIP
to retain fragments
during retrieval

Designed for **1st PASS SUCCESS** with **ALL Clot Types**

Neva™



NeVa™

Designed for **1st PASS SUCCESS**
with **ALL Clot Types**

Ideal for
**Distal M1, M2, ACA,
PCA occlusions**
Vessel diameters
2.0 – 3.5 mm

Ideal for
**MCA-M1
occlusions**
Vessel diameters
2.0 – 4.5 mm

Ideal for
**ICA tip & MCA
occlusions**
Vessel diameters
2.0 – 4.5 mm


Ideal for
**Proximal Carotid
Artery occlusions**
Vessel diameters
3.5 – 6.0 mm




NeVa M1-S
4.0 x 22 mm
30020V-MS
0.021" MC
2 Drop Zones




NeVa M1
4.0 x 30 mm
30010V-M1
0.021" MC
2 Drop Zones
Flow Restoration



NeVa T-3S
4.5 x 29 mm
VN-4529-03RR
0.021" MC
3 Drop Zones



NeVa T
4.5 x 37 mm
30040V-T
0.021" MC
2 Drop Zones
Flow Restoration



NeVa IC-3
6.0 x 44 mm
VN-6044-F3RR
0.027" MC
3 Drop Zones
Flow Restoration