

pVasc® Compatible Catheters

	pVasc 4 x 38 mm	pVasc 6 x 44 mm
Trailblazer 018 (Medtronic)	Yes	
Trailblazer 035 (Medtronic)	Yes	Yes
CXI 018 (Cook)	Yes	
CXI 035 (Cook)	Yes	Yes
Quick-Cross 018 (Phillips)	Yes	
Quick-Cross 035 (Philips)	Yes	Yes
Armada 018 Balloon (Abbott)	Yes	
Sterling 018 (Boston Scientific)	Yes	
Pacific Plus 018 PTA Catheter (Medtronic)	Yes	

Data on File.